

# HACHIKO P\*S

DE 09 58677063 Born: 08.08.2023



**Breed:** Fleckvieh  
**Breeder:** -  
**HB (DE):** 866202  
**Sire line:** HUCH

**EUROgenetik**  
**Neustadt/Aisch** freely available

<b>HASCHMICH P*S</b> DE 09 56083113 134/122/+1003-0,12-0,05	<b>HASHTAG</b> DE 09 54210676	<b>HAYABUSA</b>	<b>gGZW 134</b>
<b>TALINA</b> DE 09 54920374 123/130/+873+0,15+0,06 2/2 11238-4,86-3,80-973 2. 11792-4,99-3,79-1036	<b>GALAXI Pp*</b> DE 09 53953555	<b>MAGNET Pp*</b>	
	<b>HERKULES</b> DE 09 52599720	<b>HERZSCHLAG</b>	<b>MW 128 FW 110 FIT 107</b>
	<b>TIGI</b> DE 09 53338934 4/4 11106-4,00-3,54-838	<b>WOBLER</b>	

**BREEDING VALUES (01.04.2025)** GYB gGZW 134 (76)

GZW -1, MW +0, FW +1, FIT +0

**MILK** +981 +0,02 +43 +0,01 +36 **MW 128 (85)**

**BEEF** **FW 110 (72)**

**Net daily gain:** 107 (73) **Carcass percentage:** 114 (73) **EUROP trade class:** 102 (70)

**FITNESS** **ÖZW 131 (82)** **FIT 107 (81)**

<b>Longevity:</b> 115 (70)	<b>Udder health index:</b> 104 (83)	<b>Somatic cell count:</b> 104 (78)	<b>Mastitis:</b> 99 (60)
<b>Persistence:</b> 105 (77)	<b>Fertility index:</b> 96 (70)	<b>Claw health index:</b> 104 (71)	<b>Early fert. disorders:</b> 104 (64)
<b>Yield improvement:</b> 111 (75)	<b>Calving ease direct:</b> 103 (82)	<b>Calving ease maternal:</b> 106 (75)	<b>Cysts:</b> 95 (64)
<b>Milking speed:</b> 121 (85)	<b>Milking behaviour:</b> 100 (62)	<b>Calf vitality:</b> 104 (73)	<b>Milk fever:</b> -

**CONFORMATION** **0 daughters: 106 - 99 - 110 - 111 (82)**

Trait	76	88	100	112	124	136
<b>Frame</b> 106						
<b>Muscling</b> 99						
<b>Feet &amp; Legs</b> 110						
<b>Udder</b> 111						
Cross height 107						large
Body length 105						long
Hip width 99						wide
Body depth 104						deep
Pelvic angle 97						ascending
Hock angularity 94						straight
Hock development 101						swollen
Pasterns 101						weak
Hoof height 108						low
Fore udder length 99						short
Rear udder length 103						short
Fore udder attachment 106						loose
Central ligament 105						weak
Udder depth 108						deep
Teat length 100						short
Teat thickness 103						thin
Front teat placement 97						wide
Rear teat placement 111						wide
Rear teat direction 115						outwards
Additional teats 98						add. teats

**Free of:** ARF, B2F, DWF, F2F, F5F, TPF, ZLF  
**Genetic charac.:** P\*S, Beta casein A2A2, Kappa casein AA

